



WOMEN'S HEALTH


# NuSwab<sup>®</sup> Portfolio

A clinically relevant, cost-effective suite of tests to diagnose symptomatic patients for vaginitis/vaginosis.



# NuSwab® Portfolio


Vaginal symptoms are common in the general population and are one of the most frequent reasons for patient visits to obstetrician-gynecologists.



## NuSwab® Bacterial Vaginosis

### Result interpretation provides actionable information.

For bacterial vaginosis (BV), merely identifying the presence or absence of bacteria may not differentiate normal levels of bacteria from abnormal levels.<sup>1</sup> Numerous organisms associated with BV are also considered normal vaginal flora. Labcorp's BV assay identifies bacterial imbalance quantitatively using only three marker organisms.




## NuSwab® VG and VG+

### Targeted approach

Determining the underlying cause of vaginitis can be difficult.


- NuSwab VG combines tests for BV, *C. albicans*, *C. glabrata*, and *Trichomonas*, common causes of vaginitis.<sup>1</sup>
- NuSwab VG+ adds *Chlamydia* and *Gonorrhea* tests to NuSwab VG to aid in the identification of STI coinfection.



## NuSwab® *C. albicans* and *C. glabrata*

### Cost-effective testing with actionable results.

Two species of *Candida* comprise approximately 93% to 97%<sup>2,3</sup> of *Candida* species in vulvovaginal candidiasis according to two large US studies. Labcorp's *C. albicans* and *C. glabrata* test differentiates the two most prevalent species of *Candida*. In addition, *C. albicans* and *C. glabrata* have different – azole resistance characteristics. This differentiation may help guide appropriate therapy selection.



## NuSwab® STI

### Convenient, high-quality STI testing

NuSwab STI portfolio gives you the option to test for four common sexually transmitted infections with one vaginal swab collection.

- *Chlamydia*
- *Gonorrhea*
- *Trichomonas*
- Herpes (HSV) types 1 and 2


## Quantitative markers backed by clinical study

In a clinical study, Labcorp found excellent correlation of three BV markers (*Atopobium vaginae*, Bacterial vaginosis-associated bacterium [BVAB]-2, *Megasphaera*-1) to Nugent Gram stain score and Amsel clinical diagnosis. Labcorp found that additional BV markers studied did not improve overall test performance.<sup>4</sup>



## Organisms analyzed for concentrations and scored

Each organism's score translated into a total score.



## Result interpretation

Provider-friendly interpretation based on total score.

Total Score	Interpretation
3–6	<b>Positive</b> —indicative of bacterial vaginosis.
0–1	<b>Negative</b> —not indicative of bacterial vaginosis.
2	<b>Indeterminate</b> —unable to determine BV status. Additional clinical and laboratory data should be evaluated to establish a diagnosis.

### *Candida* species prevalence in VVC

Species	Prevalence (n=93,775) <sup>2</sup>	Prevalence (n=429) <sup>3</sup>
<i>C albicans</i>	89.0%	77.3%
<i>C glabrata</i>	7.9%	15.9%
<b>Subtotal</b>	<b>96.9%</b>	<b>93.2%</b>
<i>C parapsilosis</i>	1.7%	3.9%
<i>C krusei</i>	--	1.6%
<i>C tropicalis</i>	1.4%	1.1%
<i>C lusitaniae</i>	--	0.2%

### Target the most prevalent organisms

In two large US clinical studies, *C albicans* and *C glabrata* accounted for approximately 93% to 97%<sup>2,3</sup> of *Candida* species found. Additional *Candida* species testing is available upon request.

### Differentiate species for treatment

*C albicans* and *C glabrata* have different antifungal resistance profiles. *C glabrata* and *C krusei* have been shown to be more resistant to fluconazole. In a large study of -azole susceptibility, approximately 67% of *C glabrata* isolates demonstrated decreased susceptibility.<sup>3,5</sup>

### *Candida* species with decreased susceptibility to fluconazole<sup>3</sup>

<i>Candida</i> Species	% Resistant	% Susceptible-Dose Dependent	Total % Decreased Susceptibility	Antifungals without Decreased Susceptibility
<i>C glabrata</i>	15.2%	51.8%	67.0%	Flucytosine, Imidazoles, Nystatin
<i>C krusei</i>	Intrinsically resistant	Intrinsically resistant	Intrinsically resistant	Imidazoles, Nystatin

In the same study, *C albicans*, *C parapsilosis*, *C tropicalis*, *C lusitaniae* did not exhibit any significant decreased fluconazole susceptibility.

## Trichomonas

CDC Guidelines state that nucleic acid amplification testing (NAAT) is the recommended method to diagnose *Trichomonas vaginalis* infections. NAAT methods are high sensitive and specific and are up to 50% more sensitive for the detection of *T vaginalis* than historical methods such as wet mount microscopy, a method with poor sensitivity.<sup>6</sup>

## NuSwab® Portfolio

The NuSwab portfolio combines high-quality testing with the convenience of a single-swab collection, providing reliable and actionable information to better manage your patients. The NuSwab test menu offers a targeted approach for clinically appropriate, cost-effective care with profiles that contain fewer, select individual tests without sacrificing the content needed for comprehensive results.

	NuSwab VG	NuSwab VG+	NuSwab STI				NuSwab SELECT
Test Number	180039	180021	183160	188070	180120	180111	Individual Test List
Components	Bacterial vaginosis <i>Atopobium vaginae</i> , BVAB-2, <i>Megasphaera-1</i> <i>C albicans</i> , <i>C glabrata</i> , <i>Trichomonas</i>	Bacterial vaginosis <i>Atopobium vaginae</i> , BVAB-2, <i>Megasphaera-1</i> <i>C albicans</i> , <i>C glabrata</i> , <i>Chlamydia</i> , <i>Gonorrhea</i> , <i>Trichomonas</i>	<i>Chlamydia</i> <i>Gonorrhea</i> <i>Trichomonas</i>	<i>Chlamydia</i> <i>Gonorrhea</i> <i>Trichomonas</i> HSV 1/2	<i>Chlamydia</i> <i>Gonorrhea</i> <i>M. genitalium</i> with Reflex to Macrolide Resistance	<i>Chlamydia</i> <i>Gonorrhea</i> <i>Trichomonas</i> <i>M. genitalium</i> with Reflex to Macrolide Resistance	Bacterial vaginosis ( <b>180060</b> ) <i>C albicans</i> and <i>C glabrata</i> ( <b>180055</b> ) <i>Candida</i> Six-species Profile ( <b>180010</b> ) <i>C albicans</i> <i>C tropicalis</i> <i>C parapsilosis</i> <i>C glabrata</i> <i>C krusei</i> <i>C lusitaniae</i> <i>Chlamydia</i> / <i>Gonorrhea</i> ( <b>183194</b> ) Genital <i>Mycoplasma</i> Profile ( <b>180089</b> ) <i>M genitalium</i> , <i>M hominis</i> , <i>Ureaplasma</i> species HSV 1/2 ( <b>188056</b> ) <i>Mycoplasma genitalium</i> ( <b>180076</b> ) <i>Mycoplasma genitalium</i> , NAA, Swab With Reflex to Macrolide Resistance Testing ( <b>180092</b> ) <i>Trichomonas</i> ( <b>188052</b> )
	180042	180068					
	<b>Vaginitis (VG) With Candida (Six Species)</b> Bacterial vaginosis <i>Atopobium vaginae</i> , BVAB-2, <i>Megasphaera-1</i> <i>C albicans</i> , <i>C glabrata</i> , <i>C tropicalis</i> , <i>C parapsilosis</i> , <i>C lusitaniae</i> , <i>C krusei</i> , <i>Trichomonas</i>	<b>Vaginitis Plus (VG+) With Candida (Six Species)</b> Bacterial vaginosis <i>Atopobium vaginae</i> , BVAB-2, <i>Megasphaera-1</i> <i>C albicans</i> , <i>C glabrata</i> , <i>C tropicalis</i> , <i>C parapsilosis</i> , <i>C lusitaniae</i> , <i>C krusei</i> , <i>Chlamydia</i> , <i>Gonorrhea</i> , <i>Trichomonas</i>					
Clinical Use	Symptoms of vaginitis/vaginosis, such as discharge.	Symptoms of vaginitis/vaginosis and/or patients at risk for coinfection with Ct/Ng.	Testing patients with symptoms of multiple STIs or coinfections.				Flexibility to order any individual component.
Specimen Type	Vaginal Swab collected using the Aptima® Multitest Swab Specimen Collection Kit (preferred) or Aptima® Unisex Swab Specimen Collection Kit. Transported at room temperature.						

## Behind every test at Labcorp is a dedicated, experienced team of scientists

For the most current information regarding test options, including CPT codes, please consult the Test Menu at [womenshealth.labcorp.com/providers/sexual-health/vaginal-health](https://www.womenshealth.labcorp.com/providers/sexual-health/vaginal-health).

Labcorp's policy is to provide physicians, in each instance, with the flexibility to choose appropriate tests to assure that the convenience of ordering test combinations/profiles does not prevent physicians who wish to order a test combination/profile from making deliberate informed decisions regarding which tests are medically necessary.

All the tests offered in test combinations/profiles may be ordered individually using the Labcorp test request form.

### References

1. American College of Obstetricians and Gynecologists. Vaginitis. ACOG Practice Bulletin No. 215. *Obstet Gynecol.* Vol. 135, No. 1, January 2020.
2. Vermitsky JP, Self MJ, Chadwick SG, Trama JP, Adelson ME, Mordechai E, Gyag SE. Survey of vaginal-floral *Candida* species isolates from women of different age groups by use of species-specific PCR detection. *J Clin Microbiol.* 2008 Apr;46(4):1501-1503.
3. Richter SS, Galask R, Messer SA, Hollis RJ, Diekema DJ, Pfaller MA. Antifungal susceptibilities of *Candida* species causing vulvovaginitis and epidemiology of recurrent

cases. *J Clin Microbiol.* 2005 May;43(5):2155-2162.

4. Cartwright CP, Lembke BD, Ramachandran K, et al. Development and validation of a semiquantitative, multitarget PCR assay for diagnosis of bacterial vaginosis. *J Clin Microbiol.* 2012;50(7):2321-2329. doi:10.1128/JCM.00506-12
5. Achkar JM, Fries BC. *Candida* infections of the genitourinary tract. *Clin Micro Reviews.* 2010 Apr;23(2):253-273.
6. Centers for Disease Control and Prevention. Sexually Transmitted Diseases Treatment Guidelines, 2021. *MMWR.*

Visit the online test menu at **Labcorp.com** for additional test options and full test information, including CPT codes and specimen collection instructions.

