

Optimizing Biologic Therapy in Inflammatory Bowel Disease



DoseASSURE™, Labcorp’s portfolio of biologics monitoring assays, may help physicians maximize treatment response using a personalized, patient-specific approach

- Help aid in titrating doses or adjusting frequency to optimize effectiveness¹⁻³
- May help avoid lack of response due to under-treatment¹
- Assist in preventing and managing loss of response due to immunogenicity⁴⁻⁵
- Minimize cost to patient by avoiding unhelpful dose escalation, especially in the setting of immunogenicity¹⁻⁶

Biologic Drug Name	Primary Target*	Clinical Indications*	Test Name	Test No.
Infliximab (Remicade®; Avsola™, Renflexis®)	TNF	CD, UC **	Infliximab and Anti-Infliximab Antibody, <i>DoseASSURE™</i> IFX	503870
Adalimumab (Humira®)	TNF	CD, UC, RA	Adalimumab and Anti-Adalimumab Antibody, <i>DoseASSURE™</i> ADL	503890
Vedolizumab (Entyvio®)	α4β7 integrin	CD, UC	Vedolizumab and Anti-Vedolizumab Antibody, <i>DoseASSURE™</i> VDZ	504567
Golimumab (Simponi®)	TNF	UC, RA	Golimumab and Anti-Golimumab Antibody, <i>DoseASSURE™</i> GOL	504563
Ustekinumab (Stelara®)	IL23, IL12	CD, PA, PP	Ustekinumab and Anti-Ustekinumab Antibody, <i>DoseASSURE™</i> UST	504594
Certolizumab (Cimzia®)	TNF	CD, RA, PA, PP	Certolizumab and Anti-Certolizumab Antibody, <i>DoseASSURE™</i> CTZ	504627

*Partial listing of FDA-approved indications. TNF: tumor necrosis factor, IL: interleukin, CD: Crohn’s Disease; UC: Ulcerative Colitis, RA: Rheumatoid Arthritis, PA: Psoriatic Arthritis, PP: Plaque Psoriasis
 **Also approved for pediatric forms of CD & UC

DoseASSURE test portfolio provides tests for both drug concentration (TDM) & anti-drug antibody (immunogenicity)

Therapeutic Drug Monitoring (TDM)

- Biologics have variable pharmacokinetics.^{3,7}
- Dosing by weight and empiric dose adjustment are inefficient and suboptimal.^{3,7}
- TDM for Biologics is a valuable tool to evaluate doses and to tailor dose adjustments to your individual patient.^{3,7}
- TDM can help differentiate under-treatment from other causes of lack of response.
- Proactive dose optimization using TDM may improve clinical scores and prolong duration of anti-TNF therapy.¹
- TDM has been shown to be cost-effective and may direct more appropriate care.¹⁶

Immunogenicity Testing (Anti-drug Antibody level)

- All biologics have the potential to induce an antibody-mediated immune response.
- Close to half of IBD patients on biologic therapy may develop anti-drug antibodies.^{4,8,9}
- Anti-drug antibodies may appear as early as after the first infusion and persist for years.⁸
- Anti-drug antibodies can adversely affect the amount of drug in the body.⁸
- Sufficient drug levels (e.g. infliximab >3ug/mL), concomitant use of immunomodulating agents, and regular dosing may protect against the risk of developing anti-drug antibodies.¹⁷⁻¹⁹

Interpreting Drug Concentrations

- Detectable drug levels are associated with better clinical outcome as measured by mucosal healing, lower C-reactive protein, higher remission rate, and less relapse.^{1,2,10,11}
- Target ranges and maximally effective concentrations have not been established.³
- Optimal drug concentration depends on the desired therapeutic endpoint and may differ case by case.¹²

Drug	Normal half-life	Proposed Target Trough Concentrations§
Infliximab	7.7 to 9.5 days	3 – 7 µg/mL ¹ ; 5 -10 µg/mL ² ; >4.0 µg/mL for mucosal healing ¹² ; ≥10.0 µg/mL may be required for fistula healing ²⁰
Adalimumab	Approx 2 weeks	≥ 7.5 µg/mL ¹³ > 5.85 µg/mL ¹⁴
Vedolizumab	Approx 25 days	>30 µg/mL at week 6 ¹³ >14 µg/mL during maintenance ²⁴
Golimumab	Approx 2 weeks	≥ 4.27 µg/mL correlated with greater response and remission ²²
Ustekinumab	Approx 3 weeks	>4.5 µg/mL has been associated with greater rate of endoscopic response ²³
Certolizumab	Approx 2 weeks	≥20 µg/mL correlated to higher remission rate ¹³

§Note: These targets ranges were those used in landmark studies and do not necessarily translate into general recommendations for individual patients.

When & where to collect blood on my patients?

- The timing of sample collection is important because the drug concentration will change during the dosing interval.
- The Trough Concentration (TC) is measured at the least variable time in the dosing interval, just before the next dose (same day to within <7 days depending on the drug's normal half-life).
- During induction and maintenance phases, trough collections are usually recommended because target ranges are defined using TC.
- Blood can be drawn at any of Labcorp's nearly 2000 patient service centers located nationwide.

References

- Vande Castele N, et al. Trough Concentrations of Infliximab Guide Dosing for Patients With Inflammatory Bowel Disease. *Gastroenter*. 2015;148:1320-1329.
- Vaughn BP, et al. Proactive Therapeutic Concentration Monitoring of Infliximab May Improve Outcomes for Patients with Inflammatory Bowel Disease: Results from a Pilot Observational Study. *Inflamm Bowel Dis*. 2014;20:1996-2003.
- Vaughn BP, et al. Biologic Concentration Testing in Inflammatory Bowel Disease. *Inflamm Bowel Dis*. 2015;21:1435-1442.
- Ungar B, et al. The temporal evolution of anti-drug antibodies in patients with inflammatory bowel disease treated with infliximab. *Gut*. 2014;63:1258-1264.
- American Gastroenterological Association. Guidelines for the Identification, Assessment and Initial Medical Treatment in Crohn's Disease. <https://www.gastro.org/IBDcarepathway>
- Steenholdt C, et al. Individualised therapy is more cost-effective than dose intensification in patients with Crohn's disease who lose response to anti-TNF treatment: a randomised, controlled trial. *Gut*. 2014;63:919-927.
- Ordas, et al. Therapeutic Drug Monitoring of Tumor Necrosis Factor Antagonists in Inflammatory Bowel Disease. *Clin Gastroenterol Hepatol*. 2012;10:1079-1087.
- Steenholdt C, et al. Clinical Implications of Variations in Anti-infliximab Antibody Levels in Patients with Inflammatory Bowel Disease. *Inflamm Bowel Dis*. 2012 Vol 18(12):2209-2217
- Vande Castele N, et al. Antibody Response to Infliximab and its Impact on Pharmacokinetics can be Transient. *Am J Gastroenterol*. 2013;108:962-971.
- Maser EA, et al. Association of Trough Serum Infliximab to Clinical Outcome After Scheduled Maintenance Treatment for Crohn's Disease. *Clin Gastroenterol Hepatol*. 2006;4(10):1248-1254.
- Seovh CH, et al. Trough serum infliximab: a predictive factor of clinical outcome for infliximab treatment in acute ulcerative colitis. *Gut*. 2010;59:49-54.
- Imaeda H, et al. Relationship between serum infliximab trough levels and endoscopic activities in patients with Crohn's disease under scheduled maintenance treatment. *J Gastroenterol* 2014; 49:674-682.

Interpreting Anti-Drug Antibody Levels

- Anti-drug antibodies can impact pharmacokinetics, efficacy, and the cost effectiveness of biologics.
- Low titer antibodies may have little to no effect on drug levels or clinical outcome but evidence suggests they may lead to later development of higher titers.⁹
- In contrast, high titers of antibodies are likely to be more consequential, leading to loss of drug efficacy by preventing drug binding to TNF and/or increasing drug clearance.^{9,15}
- Anti-drug antibody positivity should be interpreted in the context of the concomitant free drug level.

Anti-Drug Antibodies	Quantitative Range	Result Interpretation
Anti-Infliximab Abs	22- 10,000+ng/mL	Antibodies are reported as Low, Intermediate or High Titer
Anti-Adalimumab Abs	25-10,000+ ng/mL	Antibodies are reported as Low, Intermediate or High Titer
Anti-Vedolizumab Abs	25-10,000+ ng/mL	Stratification into low to high titer has yet to be determined.
Anti-Golimumab Abs	20-10,000+ ng/mL	Stratification into low to high titer has yet to be determined.
Anti-Ustekinumab Abe	40-10,000+ ng/mL	Stratification into low to high titer has yet to be determined.
Anti-Certolizumab Abs	40-10,000+ ng/mL	Stratification into low to high titer has yet to be determined.

- Patient-specific clinical context must be taken into account when evaluating drug and anti-drug antibody
- Serial measurements over time may be helpful

For more information, visit labcorp.com. or call 800-444-9111.

