

Gestational Age Conversion Chart

Weeks + Day	Decimal Conversion
+1 day	0.1
+2 days	0.3
+3 days	0.4
+4 days	0.6
+5 days	0.7
+6 days	0.9

Some of Labcorp's prenatal screening assays report out gestational age in a decimal format. This guide can help you quickly convert the decimal on the report to accurately match your patient's gestational age.

Example: Gestational age of 12 weeks +3 days will report as 12.4, and 24 weeks +5 days will report as 24.7

Your laboratory partner for preconception, pregnancy and pediatric care

Maternal serum screening

Serum screening assesses risk for a range of conditions, including Down syndrome, trisomy 18 and open neural tube defects. Assays such as FirstScreen®, SequentialScreen® and Afp4® can provide a risk assessment profile specific to your patient's needs.

Preeclampsia screening Early screening saves lives.

First trimester screen – Identifies patients at higher risk of developing preeclampsia <34 weeks' gestation.

Second and third trimester test – Identifies hospitalized patients at risk of progressing to severe symptoms within 2 weeks.

For more information about our test options, visit
womenshealth.labcorp.com/providers/pregnancy

